

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 1 / 12

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

AP Density Standard Oil 30 cSt @ 40 °C

Registration number	01-2119486452-34-XXXX
IUPAC	1-Decene, homopolymer, hydrogenated
EINECS/ELINCS	500-183-1
CAS	68037-01-4

1.2 Relevant identified uses of the substance or mixture and uses advised against

1.2.1 Relevant uses

Density standard
Analytics

1.2.2 Uses advised against

None known.

1.3 Details of the supplier of the safety data sheet

Company	Anton Paar GmbH Anton-Paar-Str. 20 8054 Graz / AUSTRIA Phone +43 (0) 316 257-0 Fax +43 (0) 316 257-257 Homepage www.anton-paar.com E-mail info@anton-paar.com
---------	--

Address enquiries to

Technical information	info@anton-paar.com
Safety Data Sheet	sdb@chemiebuero.de (No dispatch of safety data sheets) Safety data sheets are available from the supplier.

1.4 Emergency telephone number

Advisory body	Call NHS 111 or a doctor
Company	+43 (0) 316 257-0

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture [REGULATION (GB) CLP]

No classification.

2.2 Label elements

The product does not require a hazard warning label in accordance with regulation CLP.

Hazard pictograms	none
Signal word	none
Hazard statements	none
Precautionary statements	none

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 2 / 12

2.3 Other hazards

Human health dangers	The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.
Environmental hazards	This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher. The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.
Other hazards	Further hazards were not determined with the current level of knowledge.

SECTION 3: Composition / Information on ingredients

3.1 Substances

The product is a substance.

Range [%]	Substance
100	1-Decene, homopolymer, hydrogenated CAS: 68037-01-4, EINECS/ELINCS: 500-183-1, Reg-No.: 01-2119486452-34-XXXX

Comment on component parts No dangerous components.
Substances of Very High Concern - SVHC: substances are not contained or are below 0.1%.

3.2 Mixtures

not applicable

SECTION 4: First aid measures

4.1 Description of first aid measures

General information	Change soaked clothing.
Inhalation	Ensure supply of fresh air. In the event of symptoms seek medical treatment.
Skin contact	When in contact with the skin, clean with soap and water. If skin irritation or rash occurs: Get medical advice/attention.
Eye contact	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
Ingestion	Do not induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious person. In the event of symptoms seek medical treatment.

4.2 Most important symptoms and effects, both acute and delayed

None known.

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Fire-fighting measures

5.1 Extinguishing media

Suitable extinguishing media	Foam, dry powder, water spray jet, carbon dioxide
Extinguishing media that must not be used	Full water jet

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 3 / 12

5.2 Special hazards arising from the substance or mixture

Carbon monoxide (CO)
Carbon dioxide (CO₂)

5.3 Advice for firefighters

Use self-contained breathing apparatus.
Collect contaminated firefighting water separately, must not be discharged into the drains.
Fire residues and contaminated firefighting water must be disposed of in accordance within the local regulations.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Use personal protective clothing.
Ensure adequate ventilation.
Remove persons to safety.
High risk of slipping due to leakage/spillage of product.

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

6.3 Methods and material for containment and cleaning up

Pick up with absorbent material (e.g. sand, universal absorbent, diatomaceous earth).
Dispose of absorbed material in accordance within the regulations.

6.4 Reference to other sections

See SECTION 8+13

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Use only in well-ventilated areas.
Avoid contact with eyes and skin. Use personal protective equipment.
The product is combustible.
Keep away from open flames, hot surfaces and sources of ignition.
Take off contaminated clothing and wash before reuse.
Wash face and/or hands before break and end of work.
Do not eat, drink or smoke when using this product.

7.2 Conditions for safe storage, including any incompatibilities

Keep only in original container.
Provide solvent-resistant and impermeable floor.
Do not store together with food and animal food/diet.
Do not store with oxidizing or self-igniting materials.
Keep container tightly closed.
Keep container in a well-ventilated place.

7.3 Specific end use(s)

See product use, SECTION 1.2

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 4 / 12

SECTION 8: Exposure controls / personal protection

8.1 Control parameters

Ingredients with occupational exposure limits to be monitored (GB)

Substance
1-Decene, homopolymer, hydrogenated
CAS: 68037-01-4, EINECS/ELINCS: 500-183-1, Reg-No.: 01-2119486452-34-XXXX
Long-term exposure: 5 mg/m ³ , OSHA PEL

Ingredients with occupational exposure limits to be monitored (EU)

not relevant

DNEL

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
There are no DNEL values established for the substance.

PNEC

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
There are no PNEC values established for the substance.

8.2 Exposure controls

Additional advice on system design	Ensure adequate ventilation on workstation.
Eye protection	If there is a risk of splashing: safety glasses (EN 166:2001)
Hand protection	The details concerned are recommendations. Please contact the glove supplier for further information. In full contact: 0.4 mm; Nitrile rubber, >480 min (EN 374-1/-2/-3).
Skin protection	Protective clothing (EN 340)
Other	Avoid contact with eyes and skin. Do not inhale aerosols.
Respiratory protection	In the event of occupational exposure limits being exceeded or of inadequate ventilation: wear appropriate respiratory protection. Short term: filter apparatus, filter A. (DIN EN 14387)
Thermal hazards	none
Delimitation and monitoring of the environmental exposition	Comply with applicable environmental regulations limiting discharge to air, water and soil.

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 5 / 12

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	liquid
Form	liquid
Color	clear colourless
Odor	odourless
Odour threshold	not applicable
pH-value	not applicable
pH-value [1%]	not applicable
Boiling point [°C]	419
Flash point [°C]	238
Flammability	not applicable
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Oxidising properties	no
Vapour pressure/gas pressure [kPa]	0.0933 (149°C)
Density [g/cm ³]	0.83 (15.6°C)
Relative density	not determined
Bulk density [kg/m ³]	not applicable
Solubility in water	insoluble
Solubility other solvents	No information available.
Partition coefficient [n-octanol/water]	> 6.5 /20°C)
Kinematic viscosity	30.5 mm ² /s (40°C)
Relative vapour density	10 (Air=1)
Evaporation speed	No information available.
Melting point [°C]	< -20
Auto-ignition temperature [°C]	354
Decomposition temperature [°C]	No information available.
Particle characteristics	not applicable

9.2 Other information

none

SECTION 10: Stability and reactivity

10.1 Reactivity

No dangerous reactions known if used as directed.

10.2 Chemical stability

Stable under normal ambient conditions (ambient temperature).

10.3 Possibility of hazardous reactions

No hazardous reactions known.

10.4 Conditions to avoid

No information available.

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 6 / 12

10.5 Incompatible materials

Oxidizing agent

10.6 Hazardous decomposition products

No dangerous reactions known if used as directed.

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 7 / 12

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute oral toxicity

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
LD50, oral, Rat, > 5000 mg/kg

Acute dermal toxicity

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
LD50, dermal, Rat, > 2000 mg/kg

Acute inhalational toxicity

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
LC50, inhalative, Rat, > 5.2mg/l 4h

Serious eye damage/irritation

Based on the available information, the classification criteria are not fulfilled.

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
Eye, Rabbit, In vivo study, non-irritating

Skin corrosion/irritation

Based on the available information, the classification criteria are not fulfilled.

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
OECD 404, no adverse effect observed (by analogy with similar substances)

Respiratory or skin sensitisation

Based on the available information, the classification criteria are not fulfilled.

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
dermal, Guinea pig, OECD 406, non-sensitizing

Specific target organ toxicity — single exposure

Based on the available information, the classification criteria are not fulfilled.

Specific target organ toxicity — repeated exposure

Based on the available information, the classification criteria are not fulfilled.

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
NOAEL, oral, Rat, 1000 mg/kg bw/day, OECD 407, no adverse effect observed

Mutagenicity

Based on the available information, the classification criteria are not fulfilled.

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
intraperitoneal, mouse, OECD 474, negativ
in vitro, OECD 471, negativ

Reproduction toxicity

Based on the available information, the classification criteria are not fulfilled.

- Fertility

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 8 / 12

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
OECD 415, no adverse effect observed (by analogy with similar substances)

- Development

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
NOAEL, oral, Rat, 1000 mg/kg bw/day, OECD 414, no adverse effect observed

Carcinogenicity Based on the available information, the classification criteria are not fulfilled.

Aspiration hazard Based on the available information, the classification criteria are not fulfilled.
On basis of test data

General remarks

The toxicity data listed pertaining to the ingredients are intended for those working in the medicinal professions, experts for occupational health and safety and toxicologists. The toxicity data pertaining to the ingredients were supplied by the manufacturers of raw materials.

11.2 Information on other hazards

Endocrine disrupting properties Contains no ingredients with endocrine-disrupting properties.

Other information none

SECTION 12: Ecological information

12.1 Toxicity

Substance
1-Decene, homopolymer, hydrogenated, CAS: 68037-01-4
LC50, (96h), Pimephales promelas, > 750 mg/l
EC50, (96h), Selenastrum capricornutum, > 1000 mg/l
EL50, (48h), Daphnia magna, > 1000 mg/l (OECD 202)

12.2 Persistence and degradability

Behaviour in environment compartments No information available.

Behaviour in sewage plant No information available.

Biological degradability This material is not expected to be readily biodegradable.
Expected to be inherently biodegradable.

12.3 Bioaccumulative potential

Accumulation in organisms is not expected.

12.4 Mobility in soil

Spillages may penetrate the soil causing ground water contamination.

12.5 Results of PBT and vPvB assessment

Based on all available information not to be classified as PBT or vPvB respectively.

12.6 Endocrine disrupting properties

Contains no ingredients with endocrine-disrupting properties.

12.7 Other adverse effects

No information available.

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 9 / 12

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Waste material must be disposed of in accordance with the Directive on waste 2008/98/EC as well as other national and local regulations. It is not possible to determine a waste code for this product in accordance with the European Waste Catalogue (EWC) since it is only possible to classify it according to how it is used by the customer. The waste code is to be determined within the EU in liaison with the waste-disposal operator.

Product

Coordinate disposal with the authorities if necessary.

Waste no. (recommended) 160509

Contaminated packaging

Uncontaminated packaging may be taken for recycling.

Waste no. (recommended) 150102
150107

SECTION 14: Transport information

14.1 UN number or ID number

Transport by land according to ADR/RID not applicable

Inland navigation (ADN) not applicable

Marine transport in accordance with IMDG not applicable

Air transport in accordance with IATA not applicable

14.2 UN proper shipping name

Transport by land according to ADR/RID NO DANGEROUS GOODS

Inland navigation (ADN) NO DANGEROUS GOODS

Marine transport in accordance with IMDG NOT CLASSIFIED AS "DANGEROUS GOODS"

Air transport in accordance with IATA NOT CLASSIFIED AS "DANGEROUS GOODS"

14.3 Transport hazard class(es)

Transport by land according to ADR/RID not applicable

Inland navigation (ADN) not applicable

Marine transport in accordance with IMDG not applicable

Air transport in accordance with IATA not applicable

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 10 / 12

14.4 Packing group

Transport by land according to ADR/RID not applicable

Inland navigation (ADN) not applicable

Marine transport in accordance with IMDG not applicable

Air transport in accordance with IATA not applicable

14.5 Environmental hazards

Transport by land according to ADR/RID no

Inland navigation (ADN) no

Marine transport in accordance with IMDG no

Air transport in accordance with IATA no

14.6 Special precautions for user

Relevant information under SECTION 6 to 8.

14.7 Maritime transport in bulk according to IMO instruments

not applicable

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

EEC-REGULATIONS 2008/98/EC 2000/532/EC; 2010/75/EU; 2004/42/EC; (EC) 648/2004; (EC) 1907/2006 (REACH); (EU) 1272/2008; 75/324/EEC ((EC) 2016/2037); (EU) 2020/878; (EU) 2016/131; (EU) 517/2014

TRANSPORT-REGULATIONS ADR (2023); IMDG-Code (2023, 41. Amdt.); IATA-DGR (2023)

NATIONAL REGULATIONS (GB): EH40/2005 Workplace exposure limits (Second edition, published December 2011); UK REACH; GB CLP.

- Observe employment restrictions for people not applicable

- VOC (2010/75/CE) 0 %

15.2 Chemical safety assessment

For this substance a chemical safety assessment has been carried out.

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 11 / 12

SECTION 16: Other information

16.1 Abbreviations and acronyms:

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route
RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses
ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure
ATE = acute toxicity estimate
CAS = Chemical Abstracts Service
CLP = Classification, Labelling and Packaging
DMEL = Derived Minimum Effect Level
DNEL = Derived No Effect Level
EC50 = Median effective concentration
ECB = European Chemicals Bureau
EEC = European Economic Community
EINECS = European Inventory of Existing Commercial Chemical Substances
EL50 = Median effective loading
ELINCS = European List of Notified Chemical Substances
EmS = Emergency Schedules
GHS = Globally Harmonized System of Classification and Labelling of Chemicals
IATA = International Air Transport Association
IBC-Code = International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
IC50 = Inhibition concentration, 50%
IMDG = International Maritime Code for Dangerous Goods
IUCLID = International Uniform Chemical Information Database
IVIS = In vitro irritation score
LC50 = Lethal concentration, 50%
LD50 = Median lethal dose
LC0 = lethal concentration, 0%
LOAEL = lowest-observed-adverse-effect level
LL50 = Median lethal loading
LQ = Limited Quantities
MARPOL = International Convention for the Prevention of Marine Pollution from Ships
NOAEL = No Observed Adverse Effect Level
NOEC = No Observed Effect Concentration
PBT = Persistent, Bioaccumulative and Toxic substance
PNEC = Predicted No-Effect Concentration
REACH = Registration, Evaluation, Authorisation and Restriction of Chemicals
STP = Sewage Treatment Plant
TLV@TWA = Threshold limit value – time-weighted average
TLV@STEL = Threshold limit value – short-time exposure limit
VOC = Volatile Organic Compounds
vPvB = very Persistent and very Bioaccumulative

16.2 Other information

Classification procedure

Anton Paar GmbH
8054 Graz

Date printed 19.04.2023, Revision 19.04.2023

Version 2.0. Supersedes version: 1.0

Page 12 / 12

Modified position

SECTION 2 been added: Further hazards were not determined with the current level of knowledge.

SECTION 2 been added: The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

SECTION 2 been added: The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

SECTION 7 been added: Provide solvent-resistant and impermeable floor.

SECTION 7 been added: Do not store with oxidizing or self-igniting materials.

SECTION 7 been added: Do not store together with food and animal food/diet.

SECTION 7 been added: Keep away from open flames, hot surfaces and sources of ignition.

SECTION 7 been added: The product is combustible.

SECTION 10 been added: Oxidizing agent

SECTION 11 been added: The toxicity data listed pertaining to the ingredients are intended for those working in the medicinal professions, experts for occupational health and safety and toxicologists. The toxicity data pertaining to the ingredients were supplied by the manufacturers of raw materials.

SECTION 11 been added: Contains no ingredients with endocrine-disrupting properties.

SECTION 11 been added: On basis of test data

SECTION 12 been added: Contains no ingredients with endocrine-disrupting properties.

SECTION 12 been added: Based on all available information not to be classified as PBT or vPvB respectively.

Copyright: Chemiebüro®

